



OptiViro[®] NK Cell Expansion Serum-free Kit Series

Product Overview

OptiVito® NK Cell Expansion Serum-free Kit Series are designed for the *in vitro* expansion of human natural killer (NK) cells, it is serum-free, xeno-free, and support high yields of functional human NK cells from peripheral blood mononuclear cells (PBMCs) and cord blood mononuclear cells (CB-MNCs).

-  **High viability**
-  **High expansion**
-  **High purity**
-  **High cytotoxicity**

NE000-N012

NK Cell Expansion Serum-free Medium P01



Basal Medium

- Serum-free, xeno-free;
- Robust NK cell expansion;
- Work with feeder or feeder-free culture system;
- GMP grade, facilitate regulatory approval.

NE000-N022

NK Cell Expansion Serum-free Kit P01



Full-cytokines Kit

- Serum-free, xeno-free;
- All-in-one, full-cytokines Kit;
- Robust NK cell expansion;
- GMP grade, facilitate regulatory approval.

NE000-N032

NK Cell Expansion Serum-free Basic Kit P01



Replenishment

- Serum-free, xeno-free;
- As a replenishment of a full-cytokines kit;
- Robust NK cell expansion;
- Work with feeder or feeder-free culture system;
- GMP grade, facilitate regulatory approval.

Product Data

1 Human NK cells from PBMCs cultured in NK Cell Expansion Serum-free Kit P01

Materials: T25 Flasks, T75 Flasks, 2L Bags

PBMCs: Seeding density at 1.3E6 cells/mL, equal to 6.5E6 cells in total.




Figure 1. High viability of NK cells.




Figure 2. High yields of total cells.




Figure 3. High yields of NK cells.

Plasma+ culture system, medium supplemented with heat-inactivated autologous plasma, obtain a density at 4.62E6 cells/mL, equal to 5.86E9 cells in total on Day 14, the culture volume is about 1.3L.

Plasma- culture system, obtain a density at 3.86E6 cells/mL, equal to 3.00E9 cells in total on Day 14, the culture volume is about 0.8L.




Figure 4. High purity of NK cells.




Figure 5. High purity of CD56+/CD16+ NK cells.

Work well without the plasma.

2 Human NK cells from CB-MNCs cultured in NK Cell Expansion Serum-free Kit P01

CB-MNCs: Seeding density at 2E6 cells/mL, equal to 2E6 cells in total.




Figure 6. High viability of NK cells.




Figure 7. High yields of total cells.




Figure 8. High purity of NK cells.

3 Human NK cells cultured demonstrate robust cytotoxicity.




Figure 9. Cytolytic function test.

Human NK cells from PBMCs, can kill target K562 cells in a dose-dependent manner.

Product Information

Product Name	Catalog No.	Amount
OptiViro® NK Cell Expansion Serum-free Medium P01	NE000-N012	1000 mL
	NE000-N011	500 mL
OptiViro® NK Cell Expansion Serum-free Kit P01	NE000-N022	1 kit (1000 mL)
	NE000-N021	1 kit (500 mL)
OptiViro® NK Cell Expansion Serum-free Basic Kit P01	NE000-N032	1 kit (1000 mL)
	NE000-N031	1 kit (500 mL)
OptiViro® NK Cell Expansion Serum-free Medium NE01 (phenol red-free)	NE000-N042	1000 mL
	NE000-N041	500 mL



Suzhou ExCell Biotechnology Co., Ltd.

☎ 400 820 5021

✉ globalsales@excellbio.com

🌐 www.excellbio.com



WeChat official account : ExCell Bio